

Short protocol

InviPrep® Fast Lysis Buffer

Isolation of DNA from Bacteria and Fungi


INVITEK
diagnostics



Zona Industrial de Tondela,
ZIM II, Lote 6
3460-070 Tondela, Portugal

techsupport@invitek.com
info@invitek.com

invitek.com



1. Shake the buffer well before use.
2. Preheat thermo shaker to 95-100°C.
3. a) **DNA purification from liquid cultures**
(e.g. culture broth)
Transfer up to 2 ml liquid sample to a Safe-Lock tube and centrifuge at max. speed (>10.000 x g) for 2 min to pellet bacterial cells. Carefully remove supernatant, do not disturb the pellet!

Add **200 µl of Fast Lysis Buffer** to the tube.

b) DNA purification from culture plates

Pick a colony from a culture plate.

Transfer cells to **200 µl of Fast Lysis Buffer** in a Safe-Lock tube.

Note: For processing difficult samples see instructions for use.

4. Resuspend by vortexing for at least 5 sec.
5. Incubate sample on thermo shaker for 8 min at 95-100°C while continuously shaking (e.g. 1400 rpm).

Note: Ensure that Safe-Lock tubes are tightly closed!

6. Cool down sample tube for 2 min at room temperature.
7. Centrifuge at max. speed (>10.000 x g) for 2 min.
8. Transfer ~50 µl to a fresh tube and avoid transferring the pellet.

Please keep the bottle of Fast Lysis Buffer together with these instructions in the outer packaging bag provided.

Note: Dilute supernatant 1:5 in DNase/RNase free water for bacterial cultures and 1:100 for fungal cultures before downstream analysis.



DOWNLOAD INSTRUCTIONS FOR USE



Read the instructions for use carefully before starting the procedure.

The instructions for use are available for download at:

[invitek.com/en/Search](https://www.invitek.com/en/Search)

invitek.com

© 2026 Invitek Diagnostics, all rights reserved.

InviPrep® Fast Lysis Buffer | 1001136100 | DE 588.03 | 2026-03-27